

# **Immunofluorescence Detection of Progesterone Receptor in Formalin-Fixed, Paraffin-Embedded Mouse Tissue**

## **Reagent and Antibody Information**

[1X Wash Buffer](#)

[3% Hydrogen Peroxide](#)

[1% BSA Diluent](#)

[1X Citrate Buffer](#)

[ProLong Gold Mounting Media](#)

### **Blocking Serum: Normal Goat Serum**

Jackson ImmunoResearch Laboratories, Inc.

West Grove, PA 19390

[www.jacksonimmuno.com](http://www.jacksonimmuno.com)

1-800-367-5296

Catalog # 005-000-121

### **Primary Antibody: Mouse Anti-Progesterone Receptor Antibody**

Beckman Coulter, Inc.

Fullerton, CA 92834

[www.beckmancoulter.com](http://www.beckmancoulter.com)

1-800-458-5060

Catalog # IM1546

### **Negative Control Serum: Purified Mouse IgG2a Isotype Control Serum**

BD Biosciences

San Jose, CA 95131

[www.bdbiosciences.com](http://www.bdbiosciences.com)

1-855-236-2772

Catalog # 550339

### **Secondary Antibody: Alexa Fluor® 488 Goat Anti-Mouse IgG (H+L) \*Highly Cross-Adsorbed\***

Life Technologies / Invitrogen

Grand Island, NY 14072

[www.invitrogen.com](http://www.invitrogen.com)

1-888-584-8929

Catalog # A-11209

## Staining Procedure

Positive Control Tissue: Female reproductive tract  
Stain localization: Nuclear

1. Deparaffinize and hydrate slides through the following solutions:

Solution	Repetitions	Time
Xylene	2 times	5 minutes
100% Ethanol	2 times	3 minutes
95% Ethanol	2 times	3 minutes
1X wash buffer	2 times	5 minutes

2. Quench endogenous peroxidase by placing the slides in 3% hydrogen peroxide for 15 minutes.
3. Rinse the slides in 2 changes of 1X wash buffer for 5 minutes each.

4. Heat-Induced Epitope Retrieval Using The Decloaker

Add 500 ml of distilled water to the pan inside the decloaker.

Place a full rack of slides into a Tissue Tek® container with 200 ml of 1X citrate buffer  
(Insert blank slides into any empty slots in the rack to ensure even heating of slides)

Place the container stably inside the pan and decloak for 5 minutes. *Maximum Pressure* \_\_\_\_\_  
Depressurize for 10 minutes.

Remove pan top and cool for 10 minutes. *Temperature Before Cooling Slides* \_\_\_\_\_

Rinse the slides in 2 changes of distilled water for 3 minutes each time.

5. Rinse the slides in 2 changes of 1X wash buffer for 5 minutes each time.
6. Block with 10% normal goat serum for 20 minutes at room temperature.  
Lot # \_\_\_\_\_ Date Reconstituted \_\_\_\_\_

**DO NOT RINSE SLIDES WITH BUFFER BEFORE ADDING PRIMARY ANTIBODY.  
ONLY WIPE EXCESS BUFFER.**

7. Apply primary antibody at a 1:150 dilution. Incubate for 1 hour at room temperature.  
Lot # \_\_\_\_\_ Exp. Date \_\_\_\_\_

For negative control slides, dilute mouse IgG2a control serum so that its IgG2a protein concentration matches that of the primary antibody (if necessary). Then make a 1:150 dilution. If the concentrations can't be matched using this method, the dilution for the negative reagent may need to be adjusted. Apply the negative and incubate for 1 hour at room temperature.

Lot # \_\_\_\_\_ Exp. Date \_\_\_\_\_

8. Rinse the slides in 2 changes of 1X wash buffer for 5 minutes each.

**From this point on, please perform the stain in the dark.**

9. Apply the Alexa 488 goat anti-mouse secondary antibody at a 1:300 dilution. Incubate for 1 hour at room temperature.

Lot # \_\_\_\_\_ Exp. Date \_\_\_\_\_

10. Rinse the slides in 2 changes of 1X Wash buffer for 5 minutes each.

11. Rinse the slides in tap water.

12. Coverslip with ProLong Gold Mounting Media (with or without DAPI).

13. Store slides at 4°C

*Updated 06/2010*