

# Detection of Synaptophysin in Formalin-Fixed, Paraffin-Embedded Rat Tissue

## **Reagent and Antibody Information**

[1X Wash Buffer](#)

[3% Hydrogen Peroxide](#)

[1% BSA Diluent](#)

[1X Citrate Buffer](#)

[Normal Rabbit IgG – Affinity Purified](#)

[DAB Chromogen](#)

[Hematoxylin](#)

### **Blocking Serum: Normal Goat Serum**

Jackson ImmunoResearch Laboratories, Inc.

West Grove, PA 19390

[www.jacksonimmuno.com](http://www.jacksonimmuno.com)

1-800-367-5296

Catalog # 005-000-121

### **Avidin / Biotin Blocking Kit**

Vector Laboratories, Inc.

Burlingame, CA 94010

[www.vectorlabs.com](http://www.vectorlabs.com)

1-800-227-6666

Catalog # SP-2001

### **Primary Antibody: Rabbit Polyclonal Anti-Synaptophysin Antibody**

Lab Vision / Thermo Fisher Scientific

Fremont, CA 94539

[www.labvision.com](http://www.labvision.com)

1-800-828-1628

Catalog # RB-1461

### **Secondary Antibody: Biotinylated Goat Anti-Rabbit IgG (H+L)**

Vector Laboratories, Inc.

Burlingame, CA 94010

[www.vectorlabs.com](http://www.vectorlabs.com)

1-800-227-6666

Catalog # BA-1000

### **Label Complex: Peroxidase-Conjugated Streptavidin SS Label**

Biogenex Laboratories

San Ramon, CA 94583

[www.biogenex.com](http://www.biogenex.com)

1-800-421-4149

Catalog # HK330-9K

## Staining Procedure

Positive Control Tissue: Pancreas (islets of Langerhans). Decreased stain intensity is observed in aged-slides (> 2 months) and over-fixed tissues.

Stain Localization: Cytoplasmic

1. Deparaffinize and hydrate slides through the following solutions:

Xylene	2 times	5 minutes
100% Ethanol	2 times	3 minutes
95% Ethanol	2 times	3 minutes
1X Wash Buffer	2 times	5 minutes

2. Quench endogenous peroxidase by placing the slides in 3% hydrogen peroxide for 15 minutes.

3. Rinse slides in 2 changes of 1X wash buffer for 5 minutes each.

4. Heat-Induced Epitope Retrieval Using The Steamer

Add distilled water to the bottom portion of the steamer.

Preheat 200ml of 1X citrate buffer in a Tissue Tek® container in the steamer between 95°C and 100°C.

Immerse a full rack of slides into the citrate buffer and place the lid back on the steamer.

(Insert blank slides into any empty slots in the rack to ensure even heating of slides)

Steam the slides for 20 minutes.

Remove container from steamer and cool for 20 minutes.

Rinse the slides in 2 changes of distilled water for 3 minutes each time.

5. Rinse the slides in 2 changes of 1X wash buffer for 5 minutes each time.

6. Block with 10% normal goat serum for 20 minutes at room temperature.

Lot # \_\_\_\_\_ Date Reconstituted \_\_\_\_\_

DO NOT RINSE SLIDES. CONTINUE TO AVIDIN-BIOTIN BLOCK.

7. Avidin / Biotin Blocking Kit

Lot # \_\_\_\_\_ Exp. Date \_\_\_\_\_ New Kit: yes / no

Apply avidin block for 15 minutes at room temperature.

Quick rinse in 1X wash buffer.

Apply biotin block for 15 minutes at room temperature.

DO NOT RINSE SECTIONS WITH BUFFER BEFORE ADDING PRIMARY ANTIBODY.  
ONLY WIPE EXCESS BLOCK.

8. Apply the primary antibody at a 1:75 dilution. Incubate for 1 hour at room temperature.

Lot # \_\_\_\_\_ Exp Date \_\_\_\_\_

For negative control slides, dilute normal rabbit IgG so that it's IgG protein concentration matches that of the primary antibody (if necessary). Then make a 1:75 dilution. If the concentrations can't be matched using this method, the dilution for the negative reagent may need to be adjusted. Apply the

negative and incubate for 1 hour at room temperature.

Lot # \_\_\_\_\_ Exp. Date \_\_\_\_\_

9. Rinse the slides in 2 changes of 1X wash buffer for 5 minutes each.

10. Apply the goat anti-rabbit secondary antibody at a 1:2000. Incubate for 15 minutes at room temperature.

Lot # \_\_\_\_\_ Date Reconstituted \_\_\_\_\_

11. Rinse the slides in 2 changes of 1X wash buffer for 5 minutes each.

12. Apply the Streptavidin SS Label. Incubate for 30 minutes at room temperature.

Lot # \_\_\_\_\_ Exp. Date \_\_\_\_\_

13. Rinse the slides in 2 changes of 1X wash buffer for 5 minutes each.

14. Apply the DAB chromogen. Incubate in the dark for 6 minutes at room temperature.

(Add 1 drop of DAB per ml of substrate)

Lot # \_\_\_\_\_ Exp. Date \_\_\_\_\_ New Kit: yes / no

15. Rinse the slides in tap water 3 minutes.

16. Counterstain with hematoxylin for 20 seconds.

17. Rinse the slides in tap water until water is clear.

18. Gently agitate slides in 1X wash buffer until the tissues turn blue.

19. Dehydrate through the following solutions:

95% Ethanol	1 time	3 minutes
100% Ethanol	3 times	3 minutes
Xylene	2 times	5 minutes

20. Coverslip

*Updated 10/28/11*