

# Detection of ED2 in Formalin-Fixed, Paraffin-Embedded Rat Tissue

## **Reagent and Antibody Information**

[1X Wash Buffer](#)

[3% Hydrogen Peroxide](#)

[1% BSA Diluent](#)

[Trypsin](#)

[DAB Chromagen](#)

[Hematoxylin](#)

### **Blocking Serum: Normal Goat Serum**

Jackson ImmunoResearch Laboratories, Inc.

West Grove, PA 19390

[www.jacksonimmuno.com](http://www.jacksonimmuno.com)

1-800-367-5296

Catalog # 005-000-121

### **Avidin / Biotin Blocking Kit**

Vector Laboratories, Inc.

Burlingame, CA 94010

[www.vectorlabs.com](http://www.vectorlabs.com)

1-800-227-6666

Catalog # SP-2001

### **Primary Antibody: Mouse Anti-Rat ED2 Antibody**

AbD Serotec, Inc.

Raleigh, NC 27604

1-919-878-7978

[www.ab-direct.com](http://www.ab-direct.com)

Catalog # MCA342R

### **Negative Control Serum: Purified Mouse IgG1 Isotype Control Serum**

BD Biosciences

San Jose, CA 95131

[www.bdbiosciences.com](http://www.bdbiosciences.com)

1-877-232-8995

Catalog # 557273

### **Secondary Antibody: Biotinylated Horse Anti-Mouse IgG (H+L)**

Vector Laboratories, Inc.

Burlingame, CA 94010

[www.vectorlabs.com](http://www.vectorlabs.com)

1-800-227-6666

Catalog # BA-2001

**Label Complex: R.T.U. Vectastain Elite ABC Reagent**

Vector Laboratories, Inc.

Burlingame, CA 94010

www.vectorlabs.com

1-800-227-6666

Catalog # PK-7100

**Staining Procedure**

Positive Control Tissue: Spleen (red pulp macrophages) and liver (Kupffer cells)

Stain Localization: Cell membrane and secreted

1. Deparaffinize and hydrate slides through the following solutions:

Solution	Repetitions	Time
Xylene	2 times	5 minutes
100% Ethanol	2 times	3 minutes
95% Ethanol	2 times	3 minutes
1X Wash Buffer	2 times	5 minutes

2. Quench endogenous peroxidase by placing slides in 3% hydrogen peroxide for 15 minutes.

3. Rinse the slides in 2 changes of 1X Wash Buffer for 5 minutes each.

4. **Proteolytic-Induced Epitope Retrieval Using Trypsin**

Incubate the slides in a 0.1% trypsin solution in a water bath at 37°C for 20 minutes.

(DO NOT add the trypsin to the 0.05M Tris-HCl • CaCl<sub>2</sub> solution until 5 minutes prior to incubation.

Trypsin loses 75% of its reactivity within 30 minutes at 37°C.)

Rinse the slides in distilled water for 1 minute to stop the enzymatic digestion.

5. Rinse the slides in 2 changes of 1X Wash Buffer for 5 minutes each.

6. Block with 10% Normal Horse Serum for 20 minutes at room temperature.

Lot # \_\_\_\_\_ Date Reconstituted \_\_\_\_\_

DO NOT RINSE SLIDES. CONTINUE TO AVIDIN-BIOTIN BLOCK.

7. **Avidin / Biotin Blocking Kit**

Lot # \_\_\_\_\_ Exp. Date \_\_\_\_\_ New Kit: yes / no

Apply avidin block for 15 minutes at room temperature.

Quick rinse in 1X Wash Buffer.

Apply biotin block for 15 minutes at room temperature.

DO NOT RINSE SECTIONS WITH BUFFER BEFORE ADDING PRIMARY ANTIBODY.  
ONLY WIPE EXCESS BLOCK.

8. Apply primary antibody at a 1:100 dilution. Incubate for 1 hour at room temperature.

Lot # \_\_\_\_\_ Date Aliquoted \_\_\_\_\_

For negative control slides, dilute the protein concentration of the mouse IgG1 serum to match that of the primary antibody, if necessary. Make a 1:100 dilution from this normalized serum, and apply to the slides. Incubate for 1 hour at room temperature.

Lot # \_\_\_\_\_ Date Reconstituted \_\_\_\_\_

9. Rinse the slides in 2 changes of 1X Wash Buffer for 5 minutes each.

10. Apply the biotinylated horse anti-mouse secondary antibody at a 1:500. Incubate for 30 minutes at room temperature.

Lot # \_\_\_\_\_ Reconstituted Date \_\_\_\_\_

11. Rinse the slides in 2 changes of 1X Wash Buffer for 5 minutes each.

12. Apply the Vectastain R.T.U Elite Label and incubate for 30 minutes at room temperature.

Exp. Date \_\_\_\_\_ New Kit: yes / no

13. Rinse slides in 2 changes of 1X Wash Buffer for 5 minutes each.

14. Apply the DAB chromagen. Incubate in the dark for 6 minutes at room temperature.

(Add 1 drop of DAB per ml of substrate)

Lot # \_\_\_\_\_ Exp. Date \_\_\_\_\_ New Kit: yes / no

15. Rinse the slides in tap water 3 minutes.

16. Counterstain with Harris Hematoxylin for 20 seconds.

17. Rinse the slides in tap water until water is clear.

18. Gently agitate slides in 1X Wash Buffer until the tissues turn blue.

19. Dehydrate through the following solutions:

Solutions	Repetitions	Time
95% Ethanol	1 time	3 minutes
100% Ethanol	3 times	3 minutes
Xylene	2 times	5 minutes

20. Coverslip