

## Overview of SNP Genotyping

Debbie Nickerson

Department of Genome Sciences  
University of Washington  
debnick@u.washington.edu


## SNP Genotyping - Overview

- Project Rationale
- Genotyping Strategies/Technical Leaps
- Data Management/Quality Control

## SNP Project Rationale

- Heritability
- Power - Number of Individuals
- Number of SNPs - Candidate Gene, Pathway, Genome  
5-10 SNPs, 400 to 1,000, 10K, 500K
- DNA requirements
- Cost

## SNP Genotyping



## SNP Typing Formats

- Microtiter Plates - Fluorescence Scale  
Low  
eg. Taqman - Good for a few markers - lots of samples - PCR prior to genotyping
- Size Analysis by Electrophoresis Medium  
eg. SNPlex - Intermediate Multiplexing reduces costs - Genotype directly on genomic DNA - new paradigm for high throughput
- Arrays - Custom or Universal High  
eg. Illumina, ParAllele, Affymetrix - Highly multiplexed - 96, 1,500 SNPs and beyond (500K+)

## Defining the scale of the genotyping project is key to selecting an approach:

	1000 individuals
5 to 10 SNPs in a candidate gene - Many approaches (expensive ~ 0.60 per SNP/genotype)	\$6,000
48 ( to 96) SNPs in a handful of candidate genes (~ 0.25 to 0.30 per SNP/genotype)	\$~29,000
384 - 1,536 SNPs - cost reductions based on scale (~0.08 - 0.15 per SNP/genotype)	\$57,600-122,880
300,000 to 500,000 SNPs defined format (~0.002 per SNP/ genotype)	\$800,000
10,000-20,000 SNPs - defined and custom formats (~0.03 per SNP/genotype)	\$>250,000


## Many Approaches to Genotype a Handful of SNPs

PCR region prior to SNP genotyping - Adds to cost  
 - Many use modified primers - the more modified, the higher the cost


- Taqman
- Single base extension - Fluorescence Polarization Sequenom - Mass Spec
- Eclipse
- Dash
- Molecular Beacons

## Taqman


Genotyping with fluorescence-based homogenous assays (single-tube assay) = 1 SNP/ tube



## Genotype Calling - Cluster Analysis



## Genotyping by Mass Spectrometry - 24



## Technological Leap - No advance PCR


Universal PCR after preparing multiple regions for analysis -





Several based on primer specific on genomic DNA followed by PCR of the ligated products - different strategies and different readouts.

SNPlex, Illumina, Parallele (Affymetrix)

Also, reduced representation - Affymetrix  
 - cut with restriction enzyme, then ligate linkers and amplify from linkers and follow by chip hybridization to read out.


## SNPlex Assay - 48 SNPs






### Arrays - High Density Genotyping Thousands of SNPs and Beyond

- "Bead" Arrays - **Illumina**
  - Manufactured by self-assembly
  - Beads identified by decoding




## Fluorescent Image of BeadArray




- ~ 3 micron diameter beads
- ~ 5 micron center-to-center
- ~50,000 features on ~1.5 mm diameter bundle
- Currently: up to 1,536 SNPs genotyped per bundle - at least 30 beads per code - many internal replicates


## Illumina Assay - 3 Primers per SNP




## Allele-Specific Extension and Ligation



## GoldenGate™ Assay Amplification





## Hybridization to Universal IllumiCode™




## BeadArray Reader

- Confocal laser scanning system
- Resolution, 0.8 micron
- Two lasers 532, 635 nm
  - Supports Cy3 & Cy5 imaging
- Sentrix Arrays (96 bundle) and Slides for 100k fixed formats







## Process Controls



## Illumina Readout for Sentrix Array > 1,000 SNPs Assayed on 96 Samples



## Multiplexed Genotyping - Universal Tag Readouts




## Parallel - Defined and Custom Formats

- Intermediate Strategy
- Multiplex ~ 20,000 SNPs
- Affymetrix readout Universal Arrays


## Parallele Technology (MIP)

### MIP Genotyping Process Overview



Molecular Inversion Probes (MIP)

### 1. Anneal



**Anneal** A mixture of Genomic DNA, up to 10,000 probes, thermostable ligase and polymerase is heat denatured and brought to annealing temperature. Two sequences located at each termini of the probe hybridize to their respective complementary sites on the genome thus forming a circular conformation with a single nucleotide gap between the termini of the probe.

### 2. Gap Fill - Polymerization



**Gap Fill polymerization** Unlabeled dATP, dCTP, dGTP or dTTP is added to each of the 4 reactions respectively. In reactions where the added nucleotide is complementary to the base being studied, DNA polymerase adds the nucleotide

### 3. Gap Fill - Ligation



**Gap Fill ligation** DNA ligase closes the gap to form a covalently closed circular molecule that encircles the genomic strand to which it is hybridized.

### 4. Exonuclease selection




**Exonuclease selection** Exonucleases are then added to digest linear probes in reactions where the added nucleotide was not complementary to the gap and excess linear probe in reactions where circular molecules were formed. The reactions are then heated to inactivate the exonucleases.

### 5. Probe release



**Probe release** The probes are then cleaved to release them from the genomic DNA



### 6. Amplification



**Amplification** The probes are amplified using common primers for all probes

Fig 3

## Affymetrix's Chip



2.5 mm  
Using Affymetrix GeneChip® Tag Array

## Whole Genome Association Strategies

Two Platforms Available Different Designs

- Affymetrix
- Illumina

## Affymetrix GeneChip Mapping 500K Array Set

## 500K: Content Optimized SNP Selection

~2,200,000 SNPs  
From Public & Perlegen

48 individuals  
Call rate, concordance


~650K SNPs

400 samples  
Call rate, accuracy  
LD

500K SNPs

- Initial Selection: 48 people
  - 2.2M SNPs
  - 25 million genotypes
  - 16 each Caucasian, African, Asian
  - All HapMap samples
- Maximize performance: Second selection over 400 people
  - 270 HapMap Samples
  - 130 diversity samples
  - Accuracy
    - HW, Mendel error, reproducibility
  - Call rates
- Maximize information content:
  - Prioritize SNPs based on LD & HapMap (Broad Institute)

## 80% genome coverage of Mapping 500K



- 500K run on 270 HapMap samples
- Pairwise  $r^2$  analysis for common SNPs ( $MAF > 0.05$ )
- Robust coverage across populations  $r^2 = 0.8$ 
  - CEPH, Asian ~66%
  - Yoruba ~45%
- 2 & 3 marker predictors (multimarker) further increase coverage


## Mapping 500K Set

- >500K SNP's
  - 2 array set
- Performance
  - 93-98% call rate range (>95% average)
  - >99.5% concordance with HapMap Genotypes, 99.9% reproducibility
- SNP lists, annotation and genotype data available without restriction at Affymetrix.com

## Illumina - Infinium I & II 10K - 300K



## Infinium II Assay Single Base Extension



## HumanHap-1 Genotyping BeadChip Content

Maximize coverage of human variation by choosing tag SNPs to uniquely identify haplotypes.

Tag SNP selection process:


1. Examine HapMap Phase I SNPs with  $MAF \geq 0.05$  in CEU
2. Bin SNPs in high LD with one another using ldSelect (Carlson, et al. 2004)
3. Select tag SNP with highest design score for each bin.



## HumanHap300 Content Strategy

- Tag SNPs
  - $r^2 \geq 0.80$  for bins containing SNPs within 10kb of genes or in evolutionarily conserved regions (ECRs)
  - $r^2 \geq 0.70$  for bins containing SNPs outside of genes or ECRs.
- Additional Content
  - ~8,000 nsSNPs
  - ~1,500 tag SNPs selected from high density SNP data in the MHC region
- Total 317,503 loci

## HumanHap300 Genomic Coverage by Population



## HumanHap300 Data Quality


127 samples  
25 trios  
15 replicates

Parameter	Percent
Call rate	99.93%
Reproducibility	>99.99%
Mendelian Inconsistencies	0.035%
Concordance with HapMap Data	99.69%

## HumanHap500 Content Strategy

- Analysis of full HapMap data set (Phase I + II) using HumanHap300 SNP list
- Fill in regions of low LD requiring higher density of tag SNPs
- Content Strategy
  - $r^2 \geq 0.80$  for bins containing SNPs within 10kb of genes or in evolutionarily conserved regions (ECRs) in CEU
  - $r^2 \geq 0.70$  for bins containing SNPs outside of genes or ECRs in CEU
  - $r^2 \geq 0.80$  for large bins ( $\geq 3$  SNPs) in CHB+JPT population
  - $r^2 \geq 0.70$  for large bins ( $\geq 5$  SNPs) in YRI population

## Preliminary HumanHap500 Genomic Coverage by Population



## Data Quality Control

- Estimating Error Rates
- Hardy Weinberg Equilibrium
- Frequency Analysis
- Missing Data

## Measuring Error Rates

- Genotype replicate samples
- Error rates generally < <1%
- Error rates are SNP specific

		Rep 1		
		CC	CT	TT
Rep 2	CC	24	1	0
	CT	0	50	0
	TT	0	0	25


## Measuring Error Rates

- Genotype replicate samples
- Absolute number of replicates is more important than percentage
  - E.g. 10% of 200?

		Rep 1		
		CC	CT	TT
Rep 2	CC	24	1	0
	CT	0	50	0
	TT	0	0	25


## Replicate samples

- Replicates can also detect sample handling errors
  - Wrong plate
  - Plate rotation



## Sample Handling Errors

- Sexing samples
- Other known genotypes
  - Blood type
  - HLA
  - Etc.



## Hardy Weinberg Equilibrium

- Given
  - $p$  = Allele 1 frequency
  - $q$  =  $1-p$
- Expectations
  - $p^2$  = frequency 11
  - $2pq$  = frequency 12
  - $q^2$  = frequency 22

## Hardy Weinberg Disequilibrium

- Heterozygote excess
  - Biologic
    - Differential survival
  - Technical
    - Nonspecific assays
    - Duplicated regions
- Homozygote excess
  - Biologic
    - Population stratification
    - Null allele
  - Technical
    - Allele dropout


## HWE Example

TABLE 2. DISTRIBUTION OF  $\alpha_2$ -ADRENERGIC-RECEPTOR VARIANTS AND  $\beta_1$ -ADRENERGIC-RECEPTOR VARIANTS AMONG CONTROLS AND PATIENTS WITH HEART FAILURE\*

ALLELES AND SUBJECTS	ALLELE FREQUENCY	P VALUE	GENOTYPE			P VALUE	ADJUSTED ODDS RATIO FOR HEART FAILURE (95% CI)†
			WT/WT	WT/Dt	Dt/Dt		
<i>α<sub>2</sub></i> D4322-325							
Black subjects							
Controls	0.411	<0.001	29/84 (34.5)	41/78 (48.8)	14/78 (16.6)	<0.001	5.65 (2.67-11.95)
Patients with heart failure	0.615		23/78 (29.5)	14/78 (17.9)	41/78 (52.6)		
White subjects							
Controls	0.038	0.01	99/105 (94.3)	4/105 (3.8)	2/105 (1.9)	0.13	3.94 (0.50-31.05)
Patients with heart failure	0.105		70/81 (86.4)	5/81 (6.2)	6/81 (7.4)		
no. total no. (%)							
Ct/Ct Cc/Cc Arg/Arg							
<i>β<sub>1</sub></i> Arg389							
Black subjects							
Controls	0.560	0.54	13/84 (15.5)	48/78 (57.1)	23/78 (27.4)	0.27	0.90 (0.44-1.84)
Patients with heart failure	0.526		19/78 (24.4)	36/78 (46.2)	23/78 (29.5)		
White subjects							
Controls	0.752	0.64	8/105 (7.6)	34/105 (32.4)	63/105 (60.0)	0.36	0.80 (0.37-1.73)
Patients with heart failure	0.741		4/81 (4.9)	34/81 (42.0)	43/81 (53.1)		

Small et al, NEJM 2002 v347 p1135

## Frequency Analysis




- Check haplotype and allele frequencies against standard populations

## Data Quality Control

- Estimate Error Rates from Replicates
- Check Hardy Weinberg Equilibrium
- Check Allele and Haplotype Frequencies
- Check Missing Data - Site specific  
Sample specific

## Detection of Outliers of the Distribution



## SNP Genotyping Summary

1. Many different genotyping approaches are available - Low to high throughput
2. Some platforms permit users to pick custom SNPs but the highest throughputs are available only in fixed contents.
3. Not all custom SNPs will work for every format. Multiple formats will be required to carry out most projects targeting specific SNPs
4. There are still trade-offs for throughput - Samples vs. SNPs
5. Costs still dictate study design.
6. Regardless of the study - Design, quality control and tracking will rule the day!  
Laboratory Information Management Systems are key in every study design  
(Key: Track - Samples,  
- Assays  
- Completion rate  
- Reproducibility/Error Analysis)