

# Detection of GFAP in Formalin-Fixed, Paraffin-Embedded Human Tissue

## **Reagent and Antibody Information**

[1X Wash Buffer](#)

[3% Hydrogen Peroxide](#)

[1% BSA Diluent](#)

[1X Citrate Buffer](#)

[DAB Chromogen](#)

[Hematoxylin](#)

**Blocking Solution: Dakocytomation Protein Block Serum-Free Ready-To-Use**

Dakocytomation Corporation

Carpinteria, CA 93013

[www.dako.com](http://www.dako.com)

1-800-235-5763

Code No. X0909

**Avidin / Biotin Blocking Kit**

Vector Laboratories, Inc.

Burlingame, CA 94010

[www.vectorlabs.com](http://www.vectorlabs.com)

1-800-227-6666

Catalog # SP-2001

**Primary Antibody: Rabbit Anti-Glial Fibrillary Acidic Protein (GFAP) Polyclonal Antibody**

Dakocytomation Corporation

Carpinteria, CA 93013

[www.dako.com](http://www.dako.com)

1-800-235-5763

Catalog # Z0334

**Negative Control Serum: Rabbit Immunoglobulin Fraction (Solid-Phase Adsorbed)**

Dakocytomation Corporation

Carpinteria, CA 93013

[www.dako.com](http://www.dako.com)

1-800-235-5763

Catalog # X0936

**Staining Kit: LSAB+ System-HRP**

Dakocytomation Corporation

Carpinteria, CA 93013

[www.dako.com](http://www.dako.com)

1-800-235-5763

Code No. K0690

**Note:** This kit includes reagents needed for the secondary antibody (link) and label complex.

## **Staining Procedure**

Positive Control Tissue: Brain – astrocytes

Stain Localization: Cytoplasmic

1. Deparaffinize and hydrate slides through the following solutions:

Solution	Repetitions	Time
Xylene	2 times	5 minutes
100% Ethanol	2 times	3 minutes
95% Ethanol	2 times	3 minutes
1X Wash Buffer	2 times	5 minutes

2. Quench endogenous peroxidase by placing the slides in 3% hydrogen peroxide for 15 minutes.

3. Rinse the slides in 2 changes of 1X wash buffer for 5 minutes each.

4. Heat-Induced Epitope Retrieval Using The Decloaker

Add 500 ml of distilled water to the pan inside the decloaker.

Place a full rack of slides into a Tissue Tek® container with 200 ml of 1X citrate buffer

(Insert blank slides into any empty slots in the rack to ensure even heating of slides)

Place the container stably inside the pan and decloak for 5 minutes. *Maximum Pressure* \_\_\_\_\_

Depressurize for 10 minutes.

Remove pan top and cool for 10 minutes. *Temperature Before Cooling Slides* \_\_\_\_\_

Rinse the slides in 2 changes of distilled water for 3 minutes each time.

5. Rinse the slides in 2 changes of 1X wash buffer for 5 minutes each.

6. Block with the Dako protein-blocking reagent for 10 minutes at room temperature.

Lot # \_\_\_\_\_ Exp Date \_\_\_\_\_

DO NOT RINSE THE SLIDES. CONTINUE TO AVIDIN-BIOTIN BLOCK.

7. Avidin / Biotin Blocking Kit

Lot # \_\_\_\_\_ Exp. Date \_\_\_\_\_ New Kit: yes / no

Apply avidin block for 15 minutes at room temperature.

Quick rinse in 1X wash buffer.

Apply biotin block for 15 minutes at room temperature.

DO NOT RINSE SLIDES WITH BUFFER BEFORE ADDING PRIMARY ANTIBODY.

ONLY WIPE EXCESS BLOCK.

8. Apply primary antibody at a 1:2000 dilution. Incubate for 30 minutes at room temperature.

Lot # \_\_\_\_\_ Exp Date \_\_\_\_\_

For negative control slides, dilute normal rabbit IgG so that it's IgG protein concentration matches that of the primary antibody (if necessary). Then make a 1:2000 dilution. If the concentrations can't be matched using this method, the dilution for the negative reagent may need to be adjusted. Apply the

negative and incubate for 30 minutes at room temperature.  
Lot # \_\_\_\_\_ Exp. Date \_\_\_\_\_

9. Rinse the slides in 2 changes of 1X wash buffer for 5 minutes each.

LSAB+ Kit  
Lot # \_\_\_\_\_ Exp Date \_\_\_\_\_

10. Apply the Link (yellow bottle) from the LSAB+ Kit. Incubate for 15 minutes at room temperature.

11. Rinse the slides in 2 changes of 1X wash buffer for 5 minutes each.

12. Apply the Label (red bottle) from the LSAB+ Kit. Incubate for 15 minutes at room temperature.

13. Rinse the slides in 2 changes of 1X wash buffer for 5 minutes each.

14. Apply the DAB chromogen. Incubate in the dark for 6 minutes at room temperature.  
(Add 1 drop of DAB per ml of substrate)  
Lot # \_\_\_\_\_ Exp Date \_\_\_\_\_ New Kit: yes / no

15. Rinse the slides in tap water 3 minutes.

16. Counterstain with hematoxylin for 20 seconds.

17. Rinse the slides in tap water until water is clear.

18. Gently agitate slides in 1X wash buffer until the tissues turn blue.

19. Dehydrate through the following solutions:

Solutions	Repetitions	Time
95% Ethanol	1 time	3 minutes
100% Ethanol	3 times	3 minutes
Xylene	2 times	5 minutes

20. Coverslip

*Updated 04/13/11*