


# UPDATE GENESPRING

Occasionally, Silicon Genetics will send you an email announcing the release of a new version of GeneSpring. Updating your computer's software will ensure that you are using the most up-to-date resources available to you. To update the GeneSpring program and documentation on your **Macintosh** or **Windows** computer, follow these steps:


1. Launch GeneSpring
2. Select **Help** from the Menu Toolbar
3. Choose the **Update GeneSpring** option




4. Be sure to place a check in the box **I agree to the above conditions.**



5. To update the online help information, click **Update Documentation**







6. The documentation will update automatically after download.
7. To update the software, click **Update Program**



8. Choose **Download** next to the platform you work with (PC, Macintosh, Solaris, or Linux)

### GeneSpring 6.2

	Platform	Installers	Instructions
	Windows	<a href="#">Download(76.5M)</a>	<a href="#">View</a>
	Mac OSX	<a href="#">Download(60.6M)</a>	<a href="#">View</a>
	Solaris	<a href="#">Download(95.9M)</a>	<a href="#">View</a>
	Linux	<a href="#">Download(98.1M)</a>	<a href="#">View</a>




For details on what's new in this release, please see our [GeneSpring 6.2 Version Notes](#).


9. Click **Yes** to install the Install Anywhere Installer




10. The installer will download



11. Select **OK** to use the English version of GeneSpring



12. Click **Next** at the **Introduction** Screen



13. Click **Next** at the **License Agreement** Screen

14. Click **Next** at the **Choose Install Folder** Screen

15. Click **Install** at the **Pre-Installation Summary** Screen

**After the installation is complete**

16. Click **Next** at the **About GeneSpring 6.2** Screen

17. Click **Done** at the **Install Complete** Screen

18. Click **Close** on the **Update Window** to complete the installation.