

Safety Data Sheet

Aflatoxins B₁, B₂, G₁, G₂, M₁, M₂, and P₁

Division of Safety
National Institutes
of Health



WARNING!

THESE COMPOUNDS MAY BE ABSORBED THROUGH THE RESPIRATORY AND INTESTINAL TRACTS. THEY MAY BE TOXIC, CARCINOGENIC, MUTAGENIC, AND TERATOGENIC. AVOID FORMATION AND BREATHING OF AEROSOLS.

LABORATORY OPERATIONS SHOULD BE CONDUCTED IN A FUME HOOD, GLOVE BOX, OR VENTILATED CABINET.

AVOID SKIN CONTACT: IF EXPOSED, WASH WITH SOAP AND WATER.

FOR EYE EXPOSURE, IRRIGATE IMMEDIATELY WITH LARGE AMOUNTS OF WATER. FOR INGESTION, DRINK WATER. INDUCE VOMITING. FOR INHALATION, REMOVE VICTIM PROMPTLY TO CLEAN AIR. ADMINISTER RESCUE BREATHING IF NECESSARY. REFER TO PHYSICIAN.

IN CASE OF LABORATORY SPILL, WEAR PROTECTIVE CLOTHING DURING CLEANUP. AVOID SKIN CONTACT OR BREATHING OF AEROSOLS. WASH DOWN AREA WITH HYPOCHLORITE SOLUTION, FOLLOWED BY 5% AQUEOUS ACETONE. SEE IARC (1980) FOR DETAILS. DISPOSE OF WASTE SOLUTIONS AND MATERIALS APPROPRIATELY.

A. Background


Aflatoxins (AF) are a class of naturally occurring compounds, produced by strains of the mold Aspergillus flavus and related species and distributed on a wide variety of foodstuffs. In pure form, they are colorless to pale yellow crystals with intense fluorescence in ultraviolet light, unstable on exposure to light and air. The term "aflatoxin" has been applied to both the natural products and some of their metabolites that retain the basic ring structure. The natural aflatoxins are divided into two classes of two compounds each on the basis of their emission,

issued 10/82


under ultraviolet light, of a blue (AFB₁, AFB₂) or green (AFG₁, AFG₂) fluorescence.* They are highly toxic in rodents, highly carcinogenic in rodents and primates, mutagenic in the Ames test, and (where data are available) teratogenic. Their use is limited to research on carcinogenicity.

B. Chemical and Physical Data


Note: The position assignments in the aflatoxin ring structure are not always consistent. The chemical names given below are based on Chemical Abstract use; in this system the double-bonded carbons of the first (left) furan ring are assigned positions 8 and 9. However, all literature reports on oxidation products (metabolites) such as the epoxide and the dihydrodiol refer to these positions as 2 and 3. To confuse the issue further, aflatoxin Q₁ is designated as 3-hydroxy-aflatoxin B₁ but in this case position 3 is consistent with Chemical Abstract nomenclature and refers to the carbon atom in meta position to the ketone group in the last (cyclopentenone) ring.




I



II



III



IV

Individual Chemical and Physical Data

Aflatoxin B₁

1. Chemical Abstract No.: 1162-65-8

Actually, this is true for AFG₁ as normally purified because of a yellow impurity. Very highly purified AFG₁ shows blue fluorescence (Lijinsky and Butler, 1966). The same may also apply to AFG₂, though there are no data.

. Synonyms:

AFB₁

Aflatoxin B

Cyclopenta(c)furo(3',2':4,5)furo(2,3-h)(1)benzopyran-
1,11-dione, 2,3,6a,9a-tetrahydro-4-methoxy- (9CI)

. Molecular

formula:

C₁₇H₁₂O₆

structure: I, R₁ = H, R₂ = CH₃

weight:

312.29

Aflatoxin B₂

Chemical Abstract No.: 7220-81-7

Synonyms:

AFB₂

Dihydroaflatoxin B₁

Cyclopenta(c)furo(3',2':4,5)furo(2,3-h)(1)benzopyran-
1,11-dione, 2,3,6a,8,9,9a-hexahydro-4-methoxy- (9CI)

Molecular

formula:

C₁₇H₁₄O₇

structure: II, R = H

weight:

314.31

Aflatoxin G₁

Chemical Abstract No.: 1165-39-5

Synonyms:

AFG₁

1H,12H-Furo(3',2':4,5)furo(2,3-h)pyrano(3,4-c)(1)benzopyran-
1,12-dione, 3,4,7a,10a-tetrahydro-5-methoxy- (9CI)

Molecular

formula:

C₁₇H₁₂O₇

structure: III, R = H

weight:

Chemical Abstract No.: 7241-98-7

Synonyms:

AFG₂

Dihydroaflatoxin G₁

1H,12H-Furo(3',2':4,5)furo(2,3h)pyrano(3,4-c)(1)benzopyran-
1,12-dione, 3,4,7a,9,10,10a-hexahydro-5-methoxy- (9CI)

Molecular

formula:

C₁₇H₁₄O₇

structure: IV, R = H

weight:

330.31

Aflatoxin M₁

Chemical Abstract No.: 6795-23-9

Synonyms

AFM₁

4-Hydroxyaflatoxin B₁

Cyclopenta(c)furo(3',2':4,5)furo(2,3-h)(1)benzopyran-
1,11-dione, 2,3,6a,9a-tetrahydro-9a-hydroxy-4-methoxy- (9CI)

Molecular

formula:

C₁₇H₁₂O₇

structure: I, R₁ = OH, R₂ = CH₃

weight:

328.29

Aflatoxin M₂

Chemical Abstract No.: 6885-57-0

Synonyms:

AFM₂

4-Hydroxyaflatoxin B₂

Cyclopenta(c)furo(3',2':4,5)furo(2,3-h)(1)benzopyran-
1,11-dione, 2,3,6a,8,9,9a-hexahydro-9a-hydroxy-4-methoxy- (9C

Molecular

formula:

C₁₇H₁₄O₇

structure: II, R = OH

weight:

Aflatoxin P₁

1. Chemical Abstract No.: 32215-02-4

2. Synonyms:

AFP₁

Cyclopenta(c)furo(3',2':4,5)furo(2,3-h)(1)benzopyran-1,11-dione, 2,3,6a,9a-tetrahydro-4-hydroxy- (9CI)

3. Molecular

formula:

C₁₆H₁₀O₆

structure: I, R₁ = H, R₂ = H

Other aflatoxins exist: for example, aflatoxin B_{2a} (Chemical Abstract No. 17878-54-5) has the same structure as aflatoxin B₂ except that an OH-group is at position 8 and aflatoxin G_{2a} (Chemical Abstract No. 20421-10-7) has the same structure as aflatoxin G₂ except that an OH-group is at position 9. Relatively little has been reported on these compounds.

Castegnaro et al. (1980) list data for the following metabolic products of aflatoxins:

1. Aflatoxin R₀ (aflatoxin F₁, aflatoxicol): reduction of cyclopentenone to the corresponding alcohol
2. Aflatoxin B_{2a} (aflatoxin B₁ hemiacetal, hydroxydihydroaflatoxin B₂)
3. Aflatoxin G_{2a} (9-hydroxyaflatoxin G₂)
4. Aflatoxin Q₁ (3-hydroxyaflatoxin B₁)

General Chemical and Physical Data

1. Density: No data.
2. Absorption spectroscopy: Naturally occurring aflatoxins and most of their identified metabolites exhibit three ultraviolet absorption maxima in the regions 214-226, 264-265, and 357-363 nm and fluorescence emission maxima at 425 or 450 nm. Individual values have been tabulated (IARC, 1976; Castegnaro et al., 1980).
3. Volatility: No data; may be considered nonvolatile.
4. Solubility: Slightly soluble in water (10-20 µg/ml); soluble in polar organic solvents such as ethanol, methanol, chloroform, dimethylsulfoxide. Aflatoxins are insoluble in nonpolar solvents, but stock and standard solutions in benzene, toluene, heptane, and cyclohexane containing a small amount of acetonitrile (Stubblefield, 1980) or n-propanol (Velasco, 1981) have been used.

5. Description, appearance: Pale yellow to colorless crystalline solids with fluorescence under UV light.

6. Boiling point: No data.

Melting point: In the range of 237-293°C with decomposition. Individual values have been tabulated (IARC, 1976). The cited values are those of the naturally occurring (optically active) aflatoxins; the synthetic (racemic) compounds have lower melting points.

7. Stability: Solid aflatoxins are stable if kept in the dark. Solutions are slowly decomposed (9-16% in 3 months) in the dark and decompose more rapidly when exposed to air and daylight or ultraviolet light; the extent varies with the solvent (Velasco, 1981). Chloroform solutions are reported to be stable for year when kept cold in the dark.

8. Chemical reactivity: The lactone ring is opened in strongly alkaline solution; the vinyl ether double bond (AFB₁, AFG₁) is reactive towards oxidizing agents (hypochlorite, perborate, ozone).

9. Flash point: No data.

10. Autoignition temperature: No data.

11. Explosive limits in air: No data.

Fire, Explosion, and Reactivity Hazard Data

1. Aflatoxins do not require special fire-fighting procedures or equipment and do not present unusual fire and explosion hazards. Because of the electrostatic nature of dry aflatoxins, fire fighters should wear full-face masks.
2. Conditions contributing to instability include exposure to heat and sunlight.
3. No incompatibilities have been reported.
4. Aflatoxins do not require nonspark equipment. When handled in flammable solvents, the precautions required for such solvents apply.

Operational Procedures

The NIH Guidelines for the Laboratory Use of Chemical Carcinogens describe operational practices to be followed when potentially carcinogenic chemicals are used in NIH laboratories. The Guidelines should be consulted to identify the proper use conditions required and specific controls to be implemented during normal and complex operations or manipulations involving AF.

1. Chemical inactivation: Validated methods have been reported (Castegnaro et al., 1980).
2. Decontamination: Turn off equipment that could be affected by AF or the materials used for cleanup. If more than 1 g has been spilled or if there is any uncertainty regarding the procedures to be followed for decontamination, call the NIH Fire Department (dial 116) for assistance. Wash surfaces with quantities of sodium hypochlorite solution. Glassware should be rinsed (in a hood) with methanol, followed by hypochlorite solution. Animal cages should be washed with hypochlorite solution. For details, see Castegnaro et al. (1980).
3. Disposal: It may be possible to decontaminate waste streams containing AF before disposal. For details, see Castegnaro et al. (1980). No waste streams containing AF shall be disposed of in sinks or general refuse. Surplus AF or chemical waste streams contaminated with AF shall be handled as hazardous chemical waste and disposed of in accordance with the NIH chemical waste disposal system. Nonchemical waste (e.g., animal carcasses and bedding) containing AF shall be handled and packaged for incineration in accordance with the NIH medical-pathological waste disposal system. Potentially infectious waste (e.g., tissue cultures) containing AF shall be disinfected by heat using a standard autoclave treatment and packaged for incineration, as above. Burnable waste (e.g., absorbent bench top liners) minimally contaminated with AF shall be handled as potentially infectious waste and packaged for incineration, as above. Absorbent materials (e.g., associated with spill cleanup) grossly contaminated shall be handled in accordance with the chemical waste disposal system. Radioactive waste containing AF shall be handled in accordance with the NIH radioactive waste disposal system.
4. Storage: Store solid materials and solutions in the cold and dark in sealed ampoules or amber screw-capped bottles or vials with Teflon cap liners. Avoid dispersal of electrostatically charged solid material while sampling.

Monitoring and Measurement Procedures Including Direct Field Measurements and Sampling for Subsequent Laboratory Analysis

1. Sampling: There are no reports on air or water sampling procedures. Various procedures have been described for sampling of agricultural commodities (Schuller et al., 1976).
2. Separation and analysis: The most commonly used procedures, as applied to food products, consist of extraction with organic solvents, separation of the extracts by column or TLC, and quantitation by densitometry or fluorescence spectroscopy. The official AOAC methods have been described (Horwitz, 1980), and solvent systems for the separation of various aflatoxins by TLC have been tabulated (Issaq and Cutchin, 1981). HPLC is the method of choice

at this time, particularly for the separation of aflatoxin metabolites from parent compounds or each other (Takahashi, 1977; Johnson et al., 1977; Colley and Neal, 1979). One application of these procedures to the analysis of aflatoxins in animal tissues has been published (Brown et al., 1973). Radioimmunoassay and enzyme-linked immunosorbent assay of aflatoxin M₁ in milk have the highest sensitivity (Pestka et al., 1981), and these methods show promise for application to other aflatoxins as well.

Biological Effects (Animal and Human)

1. Absorption: AFB₁ and AFG₁ are readily absorbed from the gastrointestinal tract in rodents and subhuman primates and from the respiratory tract in rats after intratracheal administration (Dickens et al., 1966).
2. Distribution: Oral and intraperitoneal aflatoxins or their metabolites are distributed to the liver and kidneys and probably to other organs.
3. Metabolism and excretion: Five distinct metabolic pathways for aflatoxins are known and are summarized with particular emphasis on AFB₁ (Campbell and Haynes, 1976; Colley and Neal, 1979). The first involves hydroxylation at position 4 (or position 9a in Chemical Abstracts nomenclature [refer to Section B for identification]), and the resulting products have been isolated from the milk of cows (Allcroft et al., 1966; Carnaghan et al., 1963) and urine of sheep (Holzapfel et al., 1966) fed aflatoxin-contaminated feed. Thus, AFB₁ → AFM₁ and AFB₂ → AFM₂. The second is O-demethylation (for instance, AFB₁ → AFP₁), and the resultant metabolite was isolated from the urine of monkeys dosed with AFB₁ (Dalezios et al., 1971) as the sulfate or glucuronide conjugation compound. The third pathway, so far documented only for AFB₁, is a reduction of the cyclopentanone function to "aflatoxicol," and this metabolite was identified as the major one in Sprague-Dawley rats though not in mice and monkeys, which are resistant to AFB₁-induced carcinogenesis (Wong and Hsieh, 1978). The fourth pathway consists of oxidation at the cyclopentenone ring to produce AFQ₁, which is produced in vitro by monkey and human liver microsomes and is a urinary excretion product in the monkey (Campbell and Haynes, 1976). The fifth pathway, and perhaps the most important one in terms of carcinogenicity, consists of oxidation at the 2,3 (Chem Abstract 8,9) positions of the vinyl ether double bond through the (so far only postulated) 2,3-epoxide to 2,3-dihydro-2,3-dihydroxy AFB₁. This compound has been detected as a urinary metabolite in the form of its 2-glutathionyl derivative in rats (Degen and Neumann, 1978) and is the hydrolysis product of the reaction of AFB₁ with RNA in rat and hamster liver microsomes (Swenson et al., 1974) and with DNA. Both the epoxide and aflatoxicol (Nixon et al., 1981) have been postulated to be proximate or ultimate carcinogens metabolically derived from AFB₁. Excretion is in bile, urine, and milk in unchanged form (Dann et al., 1972) or in the form of metabolites as indicated above.

4. Toxic effects: The oral LD50s of AFB₁ are 0.55, 0.62, 2.0, 2.2, 5.5, 9.0, and 10.2 mg/kg in the cat, pig, guinea pig, monkey, rat, mouse, and hamster, respectively. The only comparative toxicity studies with other aflatoxins have been carried out in day-old ducklings; oral LD50s for AFB₁, AFB₂, AFG₁, and AFG₂ are 18.2, 14.8, 39.2, and 172.5 µg per 50-g duckling (Carnaghan et al., 1963) and 16.6 and 62 µg per duckling for AFM₁ and AFM₂ (Holzapfel et al., 1966). These results indicate that metabolism by hydroxylation at the R or R₁ position does not materially alter toxicity and the less saturated compounds (structures I and III) are substantially more toxic than the others. Acute toxic effects of oral AFB₁ are on the digestive tract (emesis, anorexia); liver (jaundice, fatty degeneration, and necrosis); kidney (fatty degeneration); and central nervous system (coma, cerebral edema) in monkeys.
5. Carcinogenic effects: a. Animals. The most susceptible organ to oral AFB₁ is the liver (hepatocellular carcinomas in the rat and rhesus monkey). Other targets are kidney and colon. Inhalation of a mixture of AFB₁ and AFG₂ results in lymphatic leukemias. Feeding of several aflatoxins to pregnant rats results in hepatic carcinomas in the offspring. For a compilation of carcinogenic effects in animals, see IARC (1976).

b. Humans. It has been generally assumed from epidemiological studies in African and Asian countries that there is a positive correlation between exposure to aflatoxin-contaminated foodstuffs and incidence of hepatitis, generalized hepatotoxicity, and incidence of hepatocellular carcinoma. In recent years, the role of aflatoxins as the direct causative agent in these outbreaks of hepatocarcinomas has been questioned, since equally good correlations with exposure to hepatitis B virus can be claimed. Aflatoxins may act primarily as immunosuppressive agents, causing an increase in hepatitis B virus carriers (Lutwick, 1979; Franco et al., 1982).
6. Mutagenic and teratogenic effects: AFB₁, AFB₂, AFG₁, and AFM₁ are mutagenic in the Ames test. AFB₁ is a potent teratogen in hamsters and rats.

Emergency Treatment

1. Skin and eye exposure: For skin exposure, remove contaminated clothing and wash skin with soap and water. For eye exposure, irrigate immediately with copious quantities of running water for at least 15 minutes.
2. Ingestion: Drink plenty of water. Induce vomiting.
3. Inhalation: Remove victim promptly to clean air. Administer rescue breathing if necessary.

4. Refer to physician promptly. Undiluted bleach (5-6% NaOCl) or a mixture of sodium perborate and detergent, followed by soap and water, may be used for skin exposure. Gargles with a water solution of 1% sodium perborate and 1% sodium bicarbonate or a dentrifice containing perborate may also be used for ingestion.

References

- Allcroft, R., H. Rogers, C. Lewis, J. Nabney, and P.E. Best. 1966. Metabolism of aflatoxin in sheep: Excretion of "milk toxin." *Nature* 209:154-155.
- Brown, N.L., S. Nesheim, M.E. Stack, and G.M. Ware. 1973. Method for the determination of aflatoxin in animal tissues. *J Assoc Off Anal Chem* 56:1437-1439.
- Campbell, T.C., and J.R. Hayes. 1976. The role of aflatoxin metabolism in its toxic lesion. *Toxicol Appl Pharmacol* 35:199-222.
- Carnaghan, R.B.A., R.D. Hartley, and J. O'Kelly. 1963. Toxicity and fluorescence properties of the aflatoxins. *Nature* 200:1101.
- Castegnaro, M., D.C. Hunt, E.B. Sansone, P.L. Schuller, M.G. Siriwardana, G.M. Telling, H.P. Van Egmond, and E.A. Walker, eds. 1980. Laboratory Decontamination and Destruction of Aflatoxins B₁, B₂, G₁, G₂ in Laboratory Wastes, IARC Scientific Publications No. 37. World Health Organization, Geneva, Switzerland.
- Colley, P.J., and G.E. Neal. 1979. The analysis of aflatoxins by high-performance of liquid chromatography. *Anal Biochem* 93:409-418.
- Dalezios, J., G.N. Wogan, and S.M. Weinreb. 1971. Aflatoxin P₁: A new aflatoxin metabolite in monkeys. *Science* 171:584-585.
- Dann, R.E., L.A. Mitscher, and D. Couri. 1972. In vivo metabolism of ¹⁴C-labeled aflatoxins B₁, B₂, G₁ in rats. *Res Commun Chem Pathol Pharmacol* 3:667-675.
- Degen, G.H., and H.-G. Neumann. 1978. The major metabolite of aflatoxin B₁ in the rat is a glutathione conjugate. *Chem Biol Interact* 22:239-255.
- Dickens, F., H.E.H. Jones, and H.B. Waynforth. 1966. Oral, subcutaneous and intratracheal administration of carcinogenic lactones and related substances: The intratracheal administration of cigarette tar in the rat. *Br J Cancer* 20:134-144.
- Franco, D., D. Castaing, C. Bréchet, and J. Morin. 1982. Is aflatoxin B₁ a hepatocarcinogen in man? *Gastroenterol Clin Biol* 6:125-128.
- Holzappel, C.W., P.S. Steyn, and I.H.F. Purchase. 1966. Isolation and structure of aflatoxins M₁ and M₂. *Tetrahedron Lett* 25:2799-2803.
- Horwitz, W., ed. 1980. Official Methods of Analysis of the Association of Official Analytical Chemists, 13th edn. Chapter 26, Natural Poisons. Association of Official Analytical Chemists, Washington, DC.
- IARC, International Agency for Research on Cancer. 1976. Pages 51-72 in IARC Monographs on Evaluation of Carcinogenic Risk of Chemicals to Man, Vol. 10. World Health Organization, Geneva, Switzerland.

- Issaq, H.J., and W. Cutchin. 1981. A guide to thin-layer chromatographic systems for the separation of aflatoxins B₁, B₂, G₁ and G₂. *J Liquid Chromatogr* 4:1087-1096.
- Johnson, E., A. Abu-Shumays, and S.R. Abbott. 1977. Use of fluorescence detection in high-performance liquid chromatography. *J Chromatogr* 134:107-119.
- Lijinsky, W., and W.H. Butler. 1966. Purification and toxicity of aflatoxin G₁. *Proc Soc Exp Biol Med* 123:151-154.
- Lutwick, L.I. 1979. Relation between aflatoxin, hepatitis-B virus, and hepatocellular carcinoma. *Lancet* 1:755-757.
- Nixon, J.E., J.D. Hendricks, N.E. Pawlowski, P.M. Loveland, and R.O. Sinhuber. 1981. Carcinogenicity of aflatoxinol in Fischer 344 rats. *J Natl Cancer Inst* 66:1159-1164.
- Pestka, J.J., Y. Li, W.O. Harder, and F.S. Chu. 1981. Comparison of radioimmunoassay and enzyme-linked immunosorbent assay for determining aflatoxin M₁ in milk. *J Assoc Off Anal Chem* 64:294-301.
- Schuller, P.L., W. Horwitz, and L. Stoloff. 1976. A review of sampling plans and collaboratively studied methods of analysis for aflatoxins. *J Assoc Off Anal Chem* 59:1315-1343.
- Stubblefield, R.D. 1980. Stability and molar absorptivity of aflatoxin M₁ in acetonitrile-benzene (1+9). *J Assoc Off Anal Chem* 63:634-636.
- Swenson, D.H., J.A. Miller, and E.C. Miller. 1974. 2,3-Dihydro-2,3-dihydroxyafatoxin B₁: An acid hydrolysis product of an RNA-aflatoxin B₁ adduct formed by hamster and rat liver microsomes in vitro. *Biochem Biophys Res Commun* 53:1260-1267.
- Takahashi, D.M. 1977. High pressure liquid chromatographic determination of aflatoxins in wines and other liquid products. *J Assoc Off Anal Chem* 60:799-804.
- Velasco, J. 1981. Replacement of benzene as a solvent for aflatoxin standards. *J Am Oil Chem Soc* 58:938A-940A.
- Wong, Z.A., and D.P. Hsieh. 1978. Aflatoxinol: Major aflatoxin B₁ metabolite in rat plasma. *Science* 200:325-327.